

Technical Data Sheet

Product: Proteomic Stabilizer PROT1, 100 Test Size

Catalog Number: PROT1

Product Components:

Proteomic Stabilizer PROT1 (size: 100 Tests)

Thaw/Lyse Buffer 1, 1000X Concentrate (size: 60 milliliters)

Storage Temperature: ROOM TEMPERATURE

Shelf Life of PROT1: Recommend customer re-test every two (2) years

Shelf Life of Thaw-Lyse Buffer: 5 years as supplied

Description

Proteomic Stabilizer PROT1 is used to preserve whole blood and cell suspensions. Because applications vary we recommend that customers re-test the activity of Proteomic Stabilizer PROT1 every two (2) years to ensure there has been no loss of activity in their specific assay.

Suggested Companion Products

Thaw-Lyse Buffer, 1000X Concentrate

Lyse Buffer 2, 5X Concentrate

TruStain fcX from BioLegend

These products can be used with the following protocols. Contact Smart Tube, Inc. for copies of these protocols.

Protocol for using Bulk-Supplied Smart Tube Proteomic Stabilizer PROT1 to Bank Whole Blood Samples

Protocol Number: STP1TF-120101B

Protocol for Thawing and Lysing Erythrocytes in Whole Blood Samples Banked

with Bulk-supplied Smart Tube Proteomic Stabilizer PROT1

Protocol Number: STTLTF-130528B

Smart Tube, Inc.

Ph: (855) 397 8467

Fax: (866) 521 4037

Additional contact information: www.smarttubeinc.com/contact.htm

Any use of this product other than the permitted use without the express written authorization of Smart Tube, Inc. is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

SMART TUBE Inc.

**INNOVATIONS FOR
ON-SITE SAMPLE PREP**

WWW.SMARTTUBEINC.COM