

Technical Data Sheet

Product: Stable-Store2, 1 Liter Size

Catalog Number: STBLSTORE2-1000

Product Components:

Stable-Store V2 (size: 1 Liter)

Storage Temperature: ROOM TEMPERATURE

Shelf Life of Stable-Store: 5 years; See expiration date on label

Description

Stable-Store in combination with Stable-Lyse is used to preserve whole blood and cell suspensions.

Suggested Companion Products

STBLYSE2-250, Stable-Lyse2 250ml Size

TruStain fcX from BioLegend

These products can be used with the following protocol. Contact Smart Tube, Inc. for copies of this protocol.

Protocol for using Bulk-Supplied Stable-Lyse and Stable-Store to Bank Whole Blood Samples and Cell Suspensions. Document also contains protocol to Analyze samples previously banked with Stable-Lyse and Stable-Store Reagents.

Protocol Number: SLSSP1TF-150203

Smart Tube, Inc.

Ph: (855) 397 8467

Fax: (866) 521 4037

Additional contact information: www.smarttubeinc.com/contact.htm

Any use of this product other than the permitted use without the express written authorization of Smart Tube, Inc. is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

SMART TUBE Inc.

**INNOVATIONS FOR
ON-SITE SAMPLE PREP**

WWW.SMARTTUBEINC.COM