

Technical Data Sheet

Product: Stable-Lyse2, 250ml Size

Catalog Number: STBLYSE2-250

Product Components:

Stable-Lyse (size: 250ml)

Storage Temperature: ROOM TEMPERATURE

Shelf Life of Stable-Lyse: Recommend customer re-test every two (2) years

Description

Stable-Lyse in combination with Stable-Store is used to preserve whole blood and cell suspensions. Because applications vary we recommend that customers re-test the activity of the Stable-Lyse every two (2) years to ensure there has been no loss of activity in their specific assay.

Suggested Companion Products

STBLSTORE2-250, Stable-Store2 250ml Size

STBLSTORE2-1000, Stable-Store2 1000ml Size

TruStain fcX from BioLegend

These products can be used with the following protocol. Contact Smart Tube, Inc. for copies of this protocol.

Protocol for using Bulk-Supplied Stable-Lyse and Stable-Store to Bank Whole Blood Samples and Cell Suspensions. Document also contains protocol to Analyze samples previously banked with Stable-Lyse and Stable-Store Reagents.

Protocol Number: SLSSP1TF-150203

Smart Tube, Inc.

Ph: (855) 397 8467

Fax: (866) 521 4037

Additional contact information: www.smarttubeinc.com/contact.htm

Any use of this product other than the permitted use without the express written authorization of Smart Tube, Inc. is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

SMART TUBE Inc.

**INNOVATIONS FOR
ON-SITE SAMPLE PREP**

WWW.SMARTTUBEINC.COM